

Dialysis Cannula

Supercath™ *Clampcath*

*A gentle hemodialysis treatment
with plastic catheter tube*



- *Designed to Provide Effective Dialysis and Preserve Vascular Access*
- *Made in Japan*

CE 0197




MEDIKIT

Dialysis Cannula

Supercath™ Clampcath



Side Hole

Designed for optimal and smooth blood flow



Single Piece Jointless Plastic Catheter

Designed to prevent blood clotting and to ensure a smooth blood flow




Hub

Enables to check the blood backflow




Outer Needle

Inner Needle




Plastic Catheter Tube

Designed for smooth Catheter Tube insertion



Clamping Tube

Clamp the tube with either forceps or fingers during blood line connection



Inner Needle Tip

Designed for smooth Inner Needle insertion

Catalog Number	Catheter		Color Code	Catheter Length (mm)	No. of Side Holes	Packing
	O.D. (mm) / Gauge	I.D. (mm) / Gauge				
SP502-16 (25)	2.1 / 14G	1.7 / 16G	Pale Green	25	4	100pcs / Shelf Box
SP502-16 (33)	2.1 / 14G	1.7 / 16G		33	4	
SP502-17 (25)	1.9 / 15G	1.5 / 17G	Blue-grey	25	4	1000pcs / Shipper Box (10 Shelf Boxes)
SP502-17 (33)	1.9 / 15G	1.5 / 17G		33	4	
SP502-18 (25)	1.7 / 16G	1.3 / 18G	White	25	4	

Please read the instructions for use carefully before using the device.

*Sterilized by E.O.G (Ethylene Oxide Gas)



MEDIKIT

Medikit Co.,Ltd.

1-13-2, Yushima, Bunkyo-ku, Tokyo 113-0034, Japan
TEL:+81-3-3839-0492 | FAX:+81-3-3839-0644 | www.medikit.co.jp/english/

Manufacturer:

Togo Medikit Co., Ltd.

17148-6, Aza Kamekawa, Oaza Hichiya,
Hyuga City, Miyazaki Prefecture 883-0062, Japan

Distributed in Europe by:

Medikit Europe GmbH

Gerbermühlstrasse 7
60594 Frankfurt am Main Germany
TEL: +49-69-8700-660-60 | FAX: +49-69-8700-660-89 | www.medikit-europe.de

CE 0197

Authorized Representative in the European Community:

Bernas Medical
55 Rue Stephenson, 75018 Paris, France