

Dialysis Cannula

# Supercath™ NEO

Safety Fistula Catheter for Hemodialysis

*Fully Encased Passive  
Safety Mechanism*

*A gentle hemodialysis treatment  
with plastic catheter tube*



- *Designed to Provide Effective Dialysis and Preserve Vascular Access*
- *Made in Japan*


CE 0197



MEDIKIT

# Dialysis Cannula

# Supercath™ NEO Safety Fistula Catheter for Hemodialysis



Catalog Number	Catheter		Color Code	Catheter Length (mm)	No. of Side Holes	Packing
	O.D. (mm) / Gauge	I.D. (mm) / Gauge				
SP802-14 (25)	2.1 / 14G	1.7 / 16G	Pale Green	25	4	50pcs / Shelf Box 1000pcs / Shipper Box (20 Shelf Boxes)
SP802-14 (33)	2.1 / 14G	1.7 / 16G		33	4	
SP802-15 (25)	1.9 / 15G	1.5 / 17G	Blue-grey	25	4	
SP802-15 (33)	1.9 / 15G	1.5 / 17G		33	4	
SP802-16 (25)	1.7 / 16G	1.3 / 18G	White	25	4	

Please read the instructions for use carefully before using the device.

\*Sterilized by E.B. (Electron Beam)



**Medikit Co., Ltd.**

1-13-2, Yushima, Bunkyo-ku, Tokyo 113-0034, Japan  
TEL:+81-3-3839-0492 | FAX:+81-3-3839-0644 | www.medikit.co.jp/english/

**Manufacturer:**

**Togo Medikit Co., Ltd.**  
17148-6, Aza Kamekawa, Oaza Hichiya,  
Hyuga City, Miyazaki Prefecture 883-0062, Japan

**Distributed in Europe by:**

**Medikit Europe GmbH**  
Gerbermühlstrasse 7  
60594 Frankfurt am Main, Germany  
TEL: +49-69-8700-660-60 | FAX: +49-69-8700-660-89 | www.medikit-europe.de



**Authorized Representative in the European Community:**  
Emergo Europe B.V.  
Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands